

# Certificate of Conformity

## Pipettes

Eppendorf hereby certifies that the enclosed pipette, identified by its model and serial number, has been thoroughly inspected and validated according to the Eppendorf Quality Assurance System. This system utilizes a statistical process control for assuring the quality of design and production of each instrument. Eppendorf is an ISO 9001 certified company.

The high quality, long-term reliability of your pipette is based upon many steps of quality assurance, starting with product development and continuing through the entire production process. In addition to numerous production controls, each pipette is inspected by accredited technicians during the final test to assure that your pipette complies with its published specifications.

Accuracy of each instrument is individually checked and performance validated prior to shipment. This is indicated by a PASS status on the Check report.


A gravimetric test is part of the quality assurance system of the production of pipettes. The conditions of this gravimetric test comply with the (international) standards ISO 8655 and is also made in consideration of the standards ASTM E 1154-89 (1997).

Pipette specifications can only be guaranteed when used with original Eppendorf pipette tips. These specifications rely on the quality and consistency of the pipette system as a whole.


## Check Report

Identification:	
Model Research Series 2100 adjustable	10 - 100 µL
Serial no.	2128785
Batch	U171380L
Test Results:	
Functional Test	Pass
Volume Setting	Pass
at 10 µL Vol. setting	
Systematic error * (Inaccuracy)	+0.85 %
Random error max. (Imprecision)	0.95 %
at 100 µL Vol. setting	
Systematic error max. (Inaccuracy)	+0.14 %
Random error max. (Imprecision)	0.13 %
Inspected by:	

\*converted from the original data.



Dr. H. G. Köhn  
Managing Director



H. G. Schmidt  
Head of Quality Management

**eppendorf**  
In touch with life

Eppendorf AG · 22331 Hamburg · Germany · Phone: +49 40 538 010 · Fax: +49 40 538 01 556  
E-Mail: [eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · Internet: [www.eppendorf.com](http://www.eppendorf.com)  
Application Hotline: Phone: +49 180 366 67 89 · E-Mail: [application-hotline@eppendorf.com](mailto:application-hotline@eppendorf.com)

0012 521.833-02